

Hearing Aid Care Checklist (For Children)



Daily Care for Parents

- Wipe your child's hearing aids with a soft, dry cloth every evening
- Use a drying kit or dehumidifier box overnight (especially during monsoon or summer)
- Check for moisture or wax in the tubing / mold
- Recharge or replace batteries regularly (before bedtime or school)
- Check that both aids are working every morning before school or playtime



Weekly Maintenance

- Wash earmolds (if detachable) with warm water and mild soap
- Dry thoroughly before reattaching
- Check tubing for cracks, stiffness, or blockages
- Gently clean mic ports with a hearing aid cleaning brush



Call Your Audiologist If :

- Your child pulls the aid out frequently or seems irritated
- You hear constant whistling or feedback from the aid
- Your child doesn't respond to sounds like they normally do
- The hearing aid doesn't turn on even with a new battery
- The mold feels loose or doesn't fit snugly anymore

(Note: Children typically outgrow molds every 3–6 months)

Hearing Aid Care Checklist (For Children)



Safe Storage Tips

- Store hearing aids in a child-safe, moisture-free box every night
- Avoid keeping hearing aids in bathrooms or damp spaces
- Keep devices away from pets and toddlers
- Label the hearing aid case with your child's name (especially for school or daycare)
- Pack a spare battery and labeled case in your child's school bag



Pro Tips for Parents

- Make hearing aids part of a cheerful morning routine
- Use stickers or colored markers to help your child identify left/right aids
- Celebrate milestones – like “heard the bell” or “responded to name in noise”
- Keep a short weekly log of responses, new words, or behavioral changes

Because behind every child wearing hearing aids is a calm, consistent, and caring parent..

Soulnest is here to support every sound, every step, and every smile. 

Explore more resources, tools & 1:1 support: www.connectatsoulnest.in

